



EvoCare[®]

ECG

Compatible patient
cables with 3 and 5
leads, snap or clip
(grabber)

One-Piece Cables

3.3m cable



Leadwire Sets

0.9m cables



Trunk Cables

2.5m cable



Choosing your ECG patient cable



3 Lead Snap

Basic Trace reading
Stronger connection to
electrode press stud

[Available on request only](#)



3 Lead Clip (Grabber)

Easy connection to
electrode pad without
pressing on the patient



5 Lead Snap

Additional traces, press-on
snaps provide secure
connection reducing
noise


[Available on request only](#)



5 Lead Clip (Grabber)

Additional traces, fast
easy application to
electrode pad

AAMI One Piece ECG Cables


Models	Connector	Options	Unimed/EVO	OEM #	EVOTEC OEM
All monitor models that require the universal 6-pin round connector		3 Leads Snap	E2340S-I	852, 853, 1073/S, 9-10418-02, 8000-0025-02, 9500-0029	EA021S3I
		3 Leads Clip	E2340P-I	854, 1075/S, 8000-1005-01	EA021C3I
		5 Leads Snap	E2540S-I	854, 1075/S, 8000-1005-01	EA021S5I
		5 Leads Clip	E2540P-I	1075/P	EA021C5I

BCI (See AAMI)


Biolight One Piece ECG Cables

Models	Connector	Options	Unimed/ Evotec	OEM #	Evotec OEM
M7000, M9500 AS Modular		3 Leads Snap	E2314S-I		EA100S3I
		3 Leads Clip	E2314P-I		EA100C3I
		5 Leads Snap	E2514S-I		EA100S5I
		5 Leads Clip	E2514P-I		EA100C5I


Bionet One Piece ECG Cables

Models	Connector	Options	Unimed/ Evotec	OEM #	Evotec OEM
BM3		3 Leads Snap	E2314S-I		EA041S3I
		3 Leads Clip	E2314P-I		EA041C3I
		5 Leads Snap	E2514S-I		EA041S5I
		5 Leads Clip	E2514P-I		EA041C5I



Colin One Piece ECG Cables (Also see AAMI)

Models	Connector	Options	Unimed/ Evotec	OEM #	Evotec OEM
BP88 and BP306		3 Leads Snap	E2321S-I		EA062S3I
		3 Leads Clip	E2321P-I		EA062C3I
BP88S and Colin Press Mate Advantage		3 Leads Snap	E2320S-I	MC039-55	EA061S3I
		3 Leads Clip	E2320P-I		EA061C3I
		5 Leads Snap	E2520S-I		EA061S5I
		5 Leads Clip	E2520P-I		EA061C5I

Datascope One Piece ECG Cables

Models	Connector	Options	Unimed /Evotec	OEM #	Evotec OEM
Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Balloon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96 and System 97		3 Leads Snap	E2345S-I	0131-00-0079, 0012-00-1054-02	EA002S3I
		3 Leads Clip	E2345P-I		EA002C3I
		5 Leads Snap	E2540S-I		EA002S5I
		5 Leads Clip	E2540P-I		EA002C5I


Datascope One Piece ECG Cables

Models	Connector	Options	Unimed /Evotec	OEM #	Evotec OEM
Datascope Passport V, V12 and V21		3 Leads Snap	E2326S-I	0131-00-0079, 0012-00-1054-02	EA2326SI
		3 Leads Clip	E2326P-I		EA2326CI
		5 Leads Snap	E2526S-I		EA2526SI
		5 Leads Clip	E2526P-I		EA2526CI
Mindray Datascope Expert and Fukuda Denshi Dynascope DS-5100, DS-5100E, DS-5300W, DS-5400, DS-7001		3 Leads Snap	E2359S-I		EA068S3I
		3 Leads Clip	E2359P-I		EA068C3I
		5 Leads Snap	E2559S-I		EA068S5I
		5 Leads Clip	E2559P-I		EA068C5I


Datascope Two piece Neonatal ECG Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Datascope 2000, 2000A, 2001, 2002, 2100, 2200, 3000, 3000A, Balloon Pump, M/D3, M/D3A, M/D4, Passport, Passport IL, Passport EL, Passport XG, Passport XL, Passport 2, Passport 2LT, Passport 5L, Trio, Spectrum, Spectrum OR, System 82, System 83, System 90, System 90T, System 96 and System 97		3 lead trunk			EC002-3AI
		DIN 3 leads clip			EC024C3I

Datex Ohmeda (GE) One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
AS/3, AS/3 Compact, CS/3, S/5, S/5 Compact, Cardiocap 1, Cardiocap 2, Cardiocap 5, Light Monitor, Satellite Plus S/5, Aespire		3 Leads Snap	E2395S-I		EA010S3I
		3 Leads Clip	E2395P-I		EA010C3I
		5 Leads Snap	E2595S-I		EA010S5I
		5 Leads Clip	E2595P-I		EA010C5I

Datex Ohmeda ECG Leads

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Equipment with ASP connector for Datex Ohmeda, GE, Pro 1000, AS/3, CS/3, Cardiocap II, Cardiocap 5, Satellite Plus S/5		3 Leads Snap	EX3-90S-I	545327-HEL, 545315	EH040S3I
		3 Leads Clip	EX3-90P-I	545317-HEL, 8001958	EH040C3I
		5 Leads Snap	EX5-90S-I	545328	EH040S5I
		5 Leads Clip	EX5-90P-I	545318, 545326, 8001959-HEL	EH040C5I

Datex Ohmeda ECG Trunk Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
AS/3, CS/3, S/5, Cardiocap 1, Cardiocap 2, Cardiocap 5, Light Monitor and Satellite Plus S/5	To Monitor:  To 3 Leads Set:  To 5 Leads Set: 	3 Leads	EX-2395-I	545307, 545302	EH010-3AI
		5 Leads	EX-2595-I	545303	EH010-5AI

Datex Ohmeda DIN Leads ECG Trunk Cables

Models	Connector	Options	Part Number	OEM #	Price
DIN Safety Leadwires Trunk Cable	To Monitor:  To 3 Leads Set:  To 5 Leads Set: 	3 Leads	E-1395-I	0012-00-0991-01, 0012-00-0991-03, 1094/3, 012-108-012-0201-00	E-1395-I
		5 Leads	E-2595-I	1094/5, 012-0107-01	E-2595-I


Datex Ohmeda DIN Leads ECG Leads

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Din-Style Safety Leadwires		3 Leads Clip	E3-90P-I	0012-00-1255-09, 0012-00-1262-08, 012-0295-00, 012-0283-00, M3915A	EC024C3I
		5 Leads Clip	E5-90P-I	012-0296, 012-0284-00 012-0286-00, 012-0286-01, 012-0286-02	EC024C5I


Drager One Piece ECG Cables (Also see Siemens)

Models	Connector	Options	Unimed/Evotec	OEM #	EVOTEC OEM
UM1, UM3, UM3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, PB8800, PM8060; Series 8000, Cicero, Parameterbox		3 Leads Clip	E2377P-I		EA0573AI



Drager (Siemens) One Piece ECG Cables (Also see AAMI and Siemens)

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Sirecust Series 400, 500, 600, 700, 900, 1200, Sirecust 341, 610, 620, 630, 720, 722, 730, 732, 900, 960, 961, 1260, 1281, 1801, 400, 500 ECG Monitor, Modulo Polymed y Modulbox Sirem		3 Leads Snap	E2376S-I		
		3 Leads Clip	E2376P-I		EA052-3AI
		5 Leads Snap	E2576S-I		
		5 Leads Clip	E2576P-I		EA052-5AI


Drager (Siemens) ECG Trunk Cables (Also see AAMI and Siemens)

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Sirecust Series 400, 500, 600, 700, 900, 1200, Sirecust 341/610/620/ 630/720/722, Sirecust 730/732/900/960/961/ 1260, Sirecust 1281/ 1801/400/500, ECG Monitor, Modulo Polymed, Modulebox Sirem		3 Leads	EM-2376-I		EJ052-3AI
		5 Leads	EM-25760-I		EJ052-5AI




Drager ECG Trunk Cables (Also see Siemens)

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
For use with MultiMed Pod for Drager Infinity Vista, Vista XL, Delta, Delta XL, Gamma, Gamma XL, Gamma X XL, Kappa, Kappa XLT and Siemens 6000, 7000, 9000	To Monitor:  To 3 / 5 Leads Set: 	5 Leads, 2.5m	EM-STE5	3368391, 5950196	EQ056-5AI
MultiMed® Plus Pod for 3, 5 or 6 leads, respiration, one temperature and SpO2 for use with the following monitors: Draeger Infinity Vista, Vista XL, Delta, Delta XL, Gamma, Gamma XL, Gamma X XL, Kappa, and Kappa XL	To Monitor:  To 3 Leads Set: 	3, 5, 6 Leads	EG-STE6	MS20093	EQ091-6AI

Drager ECG Trunk Cables (Continued)


Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
SC5000, SC6000, SC7000, SC8000, SC9000	To Monitor:  To 3 Leads Set: 	3 Leads	EM-2323		EM-2323

Drager ECG Leads


Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
For use with MultiMed Pod cable for use with Drager Infinity Vista, Vista, XL, Delta, Delta XL, Gamma, Gamma XL, Gamma X XL, Kappa, Kappa XLT, Siemens 6000, 7000, 9000		3 Leads Snap	EMB3-90S-I	5956441	
		3 Leads Clip	EMB3-90P-I		
		5 Leads Snap	EM5-90S-I		
		5 Leads Clip	EM5-90P-I	5956458, 5956359	
SC5000, SC6000, SC7000, SC8000, SC9000		3 Leads Snap	EM3-90S-I		
		3 Leads Clip	EM3-90P-I		
Drager PM8060, Vitara, Parameterbox, Julian		3 Leads Clip	EMD3-90RP	M35361	
		5 Leads Clip	EMD5-90RP	M35363	
For use with Drager, MultiMed Plus, MultiMed Plus OR, and to Infinity M300 and Infinity M540 patient monitors		3 Leads Snap	EG3-90S -		
		3 Leads Clip	EG3-90P-I	MS16231	
		5 Leads Snap	EG5-90S-I		
		5 Leads Clip	EG5-90P-I	MS16546	
		6 Leads Snap	EG6-90S-I		
		6 Leads Clip	EG6-90P-I	MS16547	

Goldway (See Mindray)


Hellige One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Servomed 4, Generation, SCP 851, SCP 852, SCP 844, SCS 804, CardioServ, VICOM-sm		3 Leads Clip	E2355P-I		EA011C3I
		5 Leads Clip	E2555P-I		EA011C5I


Marquette (GE) One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Dash 1000, Dash PRO 2000, Dash 2500, Dash PRO 3000, Dash PRO 4000, Dash 5000, Eagle Systems, MAC-Lab, Solar 3000, Solar 8000, Solar 9500, Aestiva/5, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Responder 2000, Responder 3000, Corometrics 120, GEVivid 7 Ultrasound, Datex Ohmeda E Prestin modules and Unity Systems		3 Leads Snap	E2386S-I	2021141-001	EE007S3I
		3 Leads Clip	E2386P-I	2021141-001	EE007C3I
		5 Leads Snap	E2586S-I		EE007S5I
		5 Leads Clip	E2586P-I		EE007C5I

Marquette (GE) ECG Trunk Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Dash 1000, Dash Pro 2000, Dash 2500, Dash Pro MAC-Lab, Solar 3000, Solar 8000, Aestiva/5, Tram 100, Tram 200, Tram 400, Tram 500, Tram 600, Tram 800, and Carescape	To Monitor:  To 3 / 5 Leads Set:  	5 Leads	EMQ-2586-I	2017003-001, 2017003-002, 412681-005, 412931-001, E9002ZB	EE007-5I
For use with safety DIN leadwires	To Monitor:  To 3 Leads DIN Set:  To 5 Leads DIN Set: 	3 Leads	E-1386-I	CB-713006R	E-1386-I
		5 Leads	E-1586-I	CB-715006	E-1586-I



Marquette (GE) ECG Leads

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
  		3 Leads Snap	EMQ3-90S-I	416467-004, 411202-001,	EE029S3I
		3 Leads Clip	EMQ3-90P-I	412682-001, 412682-002	EE029C3I
		5 Leads Snap	EMQ5-90S-I	411202-001, 411200-001	EE029S5I
		5 Leads Clip	EMQ5-90P-I	414556-001, 414556-002, 412681-001, 412681-002, 412681-005, 414556-001,	EE029C5I
		5 Leads V2-V6 Snap	EMV5-90S-I	416468-001, 416468-003	EMV5-90S-I
		5 Leads V2-V6 Clip	EMV5-90P-I	416467-001, 416467-003	EMV5-90P-I

Medtronic (Physio Control) One Piece ECG Cables


Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Lifepak 11, Lifepak 12, Lifepak 15, Lifepak 20, Lifepak 20E	 	3 Leads Snap	E2327S-I	11110-000029, 3006218-02	EA066S3I
		3 Leads Clip	E2327P-I		EA066C3I
		5 Leads Snap	E2527S-I	11110-000066	EA066S5I
		5 Leads Clip	E2527P-I		EA066C5I

Mindray One Piece ECG Cables *(Also see Datascope)*


Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Mindray: PM9000, PM8000, PM7000 with a 6-pin connector; Goldway: UT4000A, UT4000A, UT4000Apro, UT6000A, UT4000B, UT4000C, UT4000F, UT6000F, UT4000Fpro, UT4800, Goldway G40		3 Leads Snap	E2340S-I		EA021S3I
		3 Leads Clip	E2340P-I		EA021C3I
		5 Leads Snap	E2540S-I		EA021S5I
		5 Leads Clip	E2540P-I		EA021C5I
MEC 1000, PM5000, PM6000 with a 12-pin connector, BeneView T5, T8, DPM6, DPM7 and iPM-9800. Also compatible for Penlon Intermed patient monitors		3 Leads Snap	E2312S-I	0010-30-42719	EA050S3I
		3 Leads Clip	E2312P-I		EA050C3I
		5 Leads Snap	E2512S-I		EA050S5I
		5 Leads Clip	E2512P-I		EA050C5I

Mindray – Datascope NEW – See Datascope







Mindray ECG Trunk Cables *(Also see Datascope)*

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
PM5000, PM6000, T5 and T8	To Monitor:  To 3 / 5 Leads Set:  Evotec OEM 	3 / 5 Leads 3 LEAD 5 LEAD	ER-2512-I	0010-30-42719	EJ051-3AI EJ051-5AI


Mindray ECG Leads *(Also see Datascope)*

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Mindray TMS-6016 and Tel-100 Telemetry System. Also for use with Mortara and Mennen Telemetry Systems	 Evotec OEM 	3 Leads Snap	ETM3-90S-I	115-004867-00	EJ040S3I
		3 Leads Clip	ETM3-90P-I	115-004871-00	EJ040C3I
		5 Leads Snap	ETM5-90S-I	115-004869-00	EJ040S5I
		5 Leads Clip	ETM5-90P-I	115-004873-00	EJ040C5I

Nihon Kohden One Piece ECG Cables (Also see AAMI)

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
BSM-2300 Life Scope I, BSM-2301A, BSM-2301K, BSM-2351A, BSM-4100, BSM-4102, BSM-4102A, Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, Life Scope J, BSM-9101, PVM-2701, BSM-9510 Life Scope M, BSM-9500 Life Scope S		3 Leads Snap	E2309S-I		E2309S-I
		3 Leads Clip	E2309P-I		E2309P-I
OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K, TEC-7621C, TEC-7631C, TEC-7721C, TEC-7731C, TEC-7731K, TEC-7700, Cardiolife Bifasic Defibrillator		3 Leads Snap	E2313S-I	BC-753V	EA074S3I
		3 Leads Clip	E2313P-I		EA074C3I
Merlin 78352C, Merlin 78354C, Merlin 78834C, M1001A, M1002A, Viridia, Omnicare CMA, ML65A, 66A, Heartstream and Codemaster XL		3 Leads Snap	E23029S-I		EA038S3I
		3 Leads Clip	E23029P-I		EA038C3I
		5 Leads Snap	E25029S-I		EA038S5I
		5 Leads Clip	E25029P-I		EA038C5I
OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K		3 Leads Snap	E2308S-I		EA037S3I
		3 Leads Clip	E2308P-I		EA037C3I
		5 Leads Snap	E2508S-I		EA037S5I
		5 Leads Clip	E2508P-I		EA037C5I
Life Scope OPV-1500, Life Scope N OPV-1500K		3 Leads Snap	E2316S-I		EA080S3I
		3 Leads Clip	E2316P-I		EA080C3I
PM7000, 8000, 9000		3 lead Snap			EA085S3I
		3 Lead Clip			EA085C3I
		5 Lead Snap			EA085S5I
		5 Lead Clip			EA085C3I



Nihon Kohden ECG Trunk Cables (Also see AAMI and Spacelabs)

Models	Connector	Options	Part Number	OEM #	Evotec OEM
BSM-2300 Life Scope I, BSM-2301, BSM-2351A, BSM-4100 BSM-5100 Life Scope A, Life Scope series, BSM-6301K, BSM-6501K, Life Scope J, BSM-9101, PVM-2701, Life Scope M, BSM-9800 Life Scope S, Life Scope 6700 Series	To Monitor:  To 6 Leads Set:  	6 Leads	EK-2509-I	JC-906P	EP055-6I



Nihon Kohden ECG Trunk Cables *(Also see AAMI and Spacelabs)*

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K, TEC-7100k	To Monitor:  To 3 / 5 Leads Set:  To Infant DIN Connector 	3 / 5 Leads	EK-2508-I		EC037-3I EC037-3AI
BSM-7102A, OMP-7201C, OMP-7203C, TEC-5200A (s/n<20421)	To Monitor:  To 3 / 5 Leads Set:  To DIN 	3 / 5 Leads	EK-2507-I		EC038-3I EC038-3AI




Nihon Kohden ECG Leads *(Also see AAMI and Spacelabs)*

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Cardiolife, Life Scope series 6, 8, 9, 11, 6102, 6205, 6120	  	3 Leads Snap	EKA3-90S-I		EN032S3I
		3 Leads Clip	EKA3-90P-I		EN032C3I
		5 Leads Snap	EKA5-90S-I		EN032S5I
		5 Leads Clip	EKA5-90P-I		EN032C5I
BSM-2300 Life Scope I, BSM-2301, BSM-4100 Life Scope P, BSM-5100 Life Scope A, Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, BSM-9510 Life Scope M, BSM-9800 Life Scope S, ZM-900 Series Transmitters, Life Scope 6700 Series	 	3 Leads Snap	EKB3-90S-I	BR-913P K910A, BR-913PA K910B	EP035S3I
		3 Leads Clip	EKB3-90P-I	BR-903P K911, BR-903PA K911A	EP035C3I
		6 Leads Snap	EKB6-90S-I	BR-916P K915, BR-916PA K915A	EP035S6I
		6 Leads Clip	EKB6-90P-I	BR-906P K912, BR-916PA K915A	ER035C6I


Philips One Piece ECG Cables (Also see AAMI)

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Philips HeartStart defibrillator, Philips Agilent, Philips HeartStart defibrillator, IntelliVue, IE33, 43100A, 43110A, 78354C, 78834C, 862474 C3, M1165A, M1166A, M1167A, M1204A, M1205A, M1264A, M1280A, M1722A/B, M1723A/B, M3000A, M3001A, M3046A, M2, V24E, Agilent 50XM Fetal Monitor; Solaris NT3A, NT3B, NT3C, NT3E,		3 Leads Snap	E2385S-I	M1972A, 2017003-001	EA023S3I
		3 Leads Clip	E2385P-I	M1981A	EA023C3I
		5 Leads Snap	E2585S-I	M1977A	EA023S5I
		5 Leads Clip	E2585P-I	M1986A	EA023C5I
Codemaster XL, Codemaster 100, 43100A, 43110A, 43110M, Portable, M3500B, 77000AC, Sonos/AC/AR, 77025A, 77030A, 78208A, 78210A, 78330A, 7833 1A, 78333A, 78341A, 78342A, 78345A, 78351T, 78352A/C, 78353A/B, 78534A/B/C, 78619A, 78620A, 78661A/B, 78670A, 78671A, M1723A/B, M2400A, M2406A, tems, M2425A, M2475B		3 Leads Snap	E23029S-I	M1735A	EA035S3I
		3 Leads Clip	E23029P-I		EA023C3I
		5 Leads Snap	E25029S-I	M1736A	EA023S5I
		5 Leads Clip	E25029P-I		EA023C5I



Philips ECG Trunk Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
78352C, 78354C, 78834C, Legacy HP Patient Monitors 862479 C3 Monitor 863064, 863065, 863068, SureSigns VM6, VM8, ECG M1001A/B, M1002A/B, Multi-Measurement Modules M8102A MP2 Monitor, M8105A M4735A, HeartStart XL, Defibrillator M3535A, M3536A, HeartStart MRx Defibrillator. For use with Lead Set NEW	To Monitor:  To 5 Leads Set: 	3 Leads	E-2385-I	M1669A, 989803145071	EDG073-3AI
		5 Leads	E-2585-I	M1668A, 989803145061	EDG073-5AI
Merlin, Viridia, Omnicare CMA, Heartstream, Envisor, IntelliVue, V24/V26, IE33, C Series Monitors, 78352C, 78354A/C, 78834A/B/C, M1001A, M1001B, M1002A, M1002B, M1722A, M1722B, M1723A, M1723B, M2475B, M3000A, M3001A, M4735A. For use with the Philips M1605A (LAB3-90S0) Lead Set	To Monitor:  To 3 Leads Set:  To 5 Leads Set: 	3 Leads	E-2385-I	M1500A, M1510A, M1540C, 0012- 00-0846, 4535-612-	ED023-3I
		5 Leads	E-2585-I	M1520A, M1560C, 0012-00-1002-01	ED023-5I
Codemaster XL, 782xxA Series, 783xxA Series, 77020A/AC, 77020AR, 77025A, 77040A, 78208A, 78210A, 78330A, 78332A, 78333A, 78341A, 78342A, 78345A, 78532A/B, 78534A/B/C, M1722A/B, M2400A, M2406A, M2475B, Defibrillator. For use with the Philips (LAB3-90S0) Lead Set	To Monitor:  To 3 Leads Set:  To 5 Leads Set: 	3 Leads	E-23029-I	M1580A	ED023-3I
		5 Leads	E-25029-I	M1600A	ED023-5I

Philips ECG Trunk Cables


Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
78352C, 78354C, 78834C, Legacy HP Patient Monitors 862479 C3 Monitor 863064, 863065, 863068, SureSigns VM6, VM8, ECG M1001A/B, M1002A/B, Multi- Measurement Modules M8102A MP2 Monitor, M8105A M4735A, HeartStart XL, Defibrillator M3535A, M3536A, HeartStart MRx Defibrillator. For use with Lead Set	To Monitor:  To 3 Leads Set: ORIGINAL 	3 Leads		M1669A, 989803145071	ED023C3I TO BE USED WITH ED040C3I LEADS OR ORIGINAL

Philips ECG Leads (Also see AAMI)


Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
For use with Philips M1668A ECG trunk cable INTELIVUE INTEGRATED SYSTEMS		3 Leads Snap	EP3-90S-I	M1673A	EDG040S3I
		3 Leads Clip	EP3-90P-I	M1671A, M1675A	EDG040C3I
		5 Leads Snap Limbs	EPA5-90S-I	M1644A, 989803144991	EDG040S5I
		5 Leads Clip Limbs	EPA5-90P-I	M1968A, 989803125841	
		5 Leads Snap Chest	EPB5-90S-I	M1602A	
		5 Leads Clip Chest	EPB5-90P-I	M1976A	
		3 Leads Clip	EAT3-90P-I	M2591A	
		5 Leads Snap	EAT5-90S-I	M2592A, 989803106331	
		5 Leads Clip	EAT5-90P-I	M2593A, 989803106341	
		5 Leads Clip	EHT5-90P-I	989803152051, 989803133881, 989803137251	
For use with Philips M1500A trunk cable		3 Leads Snap	EAB3-90S-I	M1605A, M1633A, M1626A, M2592A, M1613A, 0012-00-	
		3 Leads Clip	EAB3-90P-I	M1603A, M1608A, 98980310437 1,	
		5 Leads Snap	EAB5-90S-I	M1625A, 0012-00-1003-01	
		5 Leads Clip	EAB5-90P-I	M1623A, M1621A, M1629A,	

Physio Control (See AAMI and Medtronic)


Schiller One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Defigard 4000 defibrillator		3 Leads Snap	E2354S-I		E2354S-I
		3 Leads Clip	E2354P-I		E2354P-I


Siemens One Piece ECG Cables (Also see AAMI and Drager)

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Sirecust Series 400, 500, 600, 700, 900, 1200, Sirecust 341, 610, 620, 630, 720, 722, 730, 732, 900, 960, 961, 1260, 1281, 1801, 400, 500 ECG Monitor, Modulo Polymed y Modulbox Sirem		3 Leads Snap	E2376S-I		E2376S-I
		3 Leads Clip	E2376P-I		E2376P-I
		5 Leads Snap	E2576S-I		E2576S-I
		5 Leads Clip	E2576P-I		E2576P-I

Siemens ECG Trunk Cables (Also see AAMI and Drager)

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Sirecust Series 400, 500, 600, 700, 900, 1200, Sirecust 341/610/620/ 630/720/722, Sirecust 730/732/900/960/961/ 1260, Sirecust 1281/ 1801/400/500, ECG Monitor, Modulo Polymed, Modulbox Sirem		3 Leads	EM-2376-I		EJ052-3AI
		5 Leads	EM-2576-I		EJ052-5AI

Welch Allyn One Piece ECG Cables

Models	Connector	Options	Unimed/Evotec	OEM #	Evotec OEM
Propaq 100 Series, Propaq 200 Series, Propaq Encore, Propaq CS, Ultra Smart, Atlas, PIC 30, PIC 40 and PIC 50		3 Leads Snap	E2340S-I		EA021S-3I
		3 Leads Clip	E2340P-I		ES021C-3I
		5 Leads Snap	E2540S-I		ES021S-5I
		5 Leads Clip	E2540P-I		ES021C-5I

Zoll (See AAMI)