



EvoCare[®]

IBP

IBP Adapters

3m cable, connects to disposable transducer



Compatible IBP adapters, converter cables and disposable transducers

IBP Converters

0.3m cable



Transducers

Sterile, single-use, disposable




Pressure Infusion Bags

Disposable and reusable



Choosing your IBP adapter cable


IBP adapters are available with the following connectors for use with all of the major disposable transducer brands:




Argon




B. Brown




BD



Edwards



Abbott




Utah



Medex Logical

6-Pin IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
IBP adapter cable for Bard H8000 Intra Aortic Balloon Pump (No Cal); BCI 9100 Series and 9200 Advisor Series; Adecon, Critikon; Digicare; Dixtal DX-2020, DX-2021, and DX-2023; Genera 710C; Hewlett Packard 5microvolt, (Option 002) 8805, 78205; Invivo Corporation/Medical Data Electronics (NDE), Inc. TVS Omnitrak, Omnitrak 4000, Millenia, Model 3550 Escort M12, Escort 100 - 200 - and 300 Series, MDE, Series 3100 - 3500; Ivy Biomedical 700, 400 Series; Medical Research Laboratories - Welch Allyn Portable Invasive Care (PIC); Mindray / Datascope Passport; Mennen; Nellcor Escort; Omni-Trak 4000; PaceTech VitalMax 4000 Modular; Physio Control Lifepak 12; Smiths Medical; Spacelabs 300, 400, 500, 600, Alpha 7, Alpha PC, 90308, 90311, 90331, 90402, 90403, 90470, 90496*, 90603A, 91496*, PC Express, 700 Surgical Series, MOD-CD, Tel/Vitatek, Scout TPC *90496 and 91496 require Spacelabs 700-0028-00 rectangular jumper cable; Tecktronix 400 Series; Vitalert 2000; Viatek; Welch Allyn Propaq 102, 104, 106 Acuity, 106EL, Encore; Zoll Defibrillators		Argon	EC-6P-AG		EV6P-AR
		B. Braun	EC-6P-BB		EV6P-BB
		BD	EC-6P-BD	001C-30-70758, 684085	EV6P-BD
		Edwards	EC-6P-ED	0010-21-43094	EV6P-ED
		Medex Abbott	EC-6P-MX		EV6P-MA
		Utah	EC-6P-UT	650-208	EV6P-UT
		Medex Logical	EC-6P-MX1		EV6P-ML

ADECON (See 6-Pin IBP Adapter)

AEONMED (See 6-Pin IBP Adapter)

Air Shields (See Datex Ohmeda)


BCI (See 6-Pin IBP Adapter)

Bionics (See 6-Pin IBP Adapter)

Burdick (See Datex Ohmeda)

Colin (See Datascope)


Biolight BLT IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Anyview range V6, A2, A3, A5, A6		Argon	EC-BLT-AG		EVBLT-AR
		B. Braun	EC-BLT-BB		EVBLT-BB
		BD	EC-BLT-BD		EVBLT-BD
		Edwards	EC-BLT-ED		EVBLT-ED
		Medex Abbott	EC-BLT-MX		EVCBLT-MA
		Utah	EC-BLT-UT		EVBLT-UT
		Medex Logical	EC-BLT-MX1		EVBLT-ML

Bionics IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Bpm1200		Argon	EC-BLT-AG		EVBLT-AR
		B. Braun	EC-BLT-BB		EVBLT-BB
		BD	EC-BLT-BD		EVBLT-BD
		Edwards	EC-BLT-ED		EVBLT-ED
		Medex Abbott	EC-BLT-MX		EVCBLT-MA
		Utah	EC-BLT-UT		EVBLT-UT
		Medex Logical	EC-BLT-MX1		EVBLT-ML

Criticare IBP Adapter



Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Criticare Matrix, Scholar 2200		Argon	EC-CSI-AG		EVCSI-AR
		B. Braun	EC-CSI-BB		EVCSI-BB
		BD	EC-CSI-BD	684175	EVCSI-BD
		Edwards	EC-CSI-ED		EVCSI-ED
		Medex Abbott	EC-CSI-MX		EVCSI-MA
		Utah	EC-CSI-UT		EVCSI-UT
		Medex Logical	EC-CSI-MX1		EVCSI-ML

Critikon (See 6-Pin IBP Adapter)


Colin (See Datascope)

Contec (See 6-Pin IBP Adapter)


Comen IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
C100, 90170, 60, 50, NC8		Argon			EVCO-AR
		B. Braun			EVCO-BB
		BD			EVCO-BD
		Edwards			EVCO-ED
		Medex Abbott			EVCO-MA
		Utah			EVCO-UT
		Medex Logical			EVCO-ML
C80		Argon			EVCO1-AR
		B. Braun			EVCO1-BB
		BD			EVCO1-BD
		Edwards			EVCO1-ED
		Medex Abbott			EVCO1-MA
		Utah			EVCO1-UT
		Medex Logical			EVCO1-ML


Datascope IBP Adapter (Also see 6-Pin IBP Adapter)

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
American Optical 1131 Module; Colin 408BP, 400 Series; CWE Research Lab Monitor; Datascope System IA Balloon Pump 95, 97, 98, 98XT, CS100, CS300, 2000, 2000A, 2001, 2002, 3000, Passport, Passport II, Passport XG, Spectrum, Spectrum OR, Trio; Hewlett Packard Sandborn, INEX; Kontron (Roche) 2106, 7320, 7324 (BP Module for 7046 & 7048); Marquette 3100 Series		Argon	EC-DT-AG		EVDT-AR
		B. Braun	EC-DT-BB		EVDT-BB
		BD	EC-DT-BD	684078	EVDT-BD
		Edwards	EC-DT-ED		EVDT-ED
		Medex Abbott	EC-DT-MX		EVDT-MA
		Utah	EC-DT-UT	650-204	EVDT-UT
		Medex Logical	EC-DT-MX1		EVDT-ML



Datascope IBP Adapter (Continued)

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Datascope Passport V, V12, V21		Argon	EC-DT1-AG		EVDT1-AR
		B. Braun	EC-DT1-BB		EVDT1-BB
		BD	EC-DT1-BD	684078	EVDT1-BD
		Edwards	EC-DT1-ED		EVDT1-ED
		Medex Abbott	EC-DT1-MX		EVDT1-MA
		Utah	EC-DT1-UT	650-204	EVDT1-UT
		Medex Logical	EC-DT1-MX1		EVDT1-ML


Datex Ohmeda (GE) IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Air-Shields AS-100 (141 Module), AS-441; Burdick M500 Series M544 Module (M560/M565/M575 Mainframe); Diascope Traveler/North American; GE Datex Ohmeda AS/3, AS/5, ESTP Module, CC-104, Cardiocap CM-104, Cardiocap I and II, Cardiocap 5; Kone M500 Series M544 Module (M560/ M565/ M575 Mainframe); Puritan - Bennett PB 240, PB 245, PB 245P; Simonson & Weil AS100 Series, Athena System, 8000 Quadrascope		Argon	EC-DX-AG		EVDA-AR
		B. Braun	EC-DX-BB		EVDA-BB
		BD	EC-DX-BD		EVDA-BD
		Edwards	EC-DX-ED	684104	EVDA-ED
		Medex Abbott	EC-DX-MX		EVDA-MA
		Utah	EC-DX-UT	650-217	EVDA-UT
		Medex Logical	EC-DX1-MX1		EVDA-ML

Drager IBP Adapter **Also see converter cables below**


Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Siemens Axiom Sensis XP with Hemopod, 3rds Generation Sirecust 400 Series, 630, 700, 732, 960, 961, 1241, 1260, 1280, 1281, 4041, SC6000, SC6002, SC6002 XL, SC6802 XL, Gamma XL, Vista XL, SC7000, SC8000, SC9000, SC9000 XL, Delta, Delta XL		Argon	EC-SM1-AG		EVDR-AR
		B. Braun	EC-SM1-BB		EVDR-BB
		BD	EC-SM1-BD	684082	EVDR-BD
		Edwards	EC-SM1-ED		EVDR-ED
		Medex Abbott	EC-SM1-MX		EVDR-MA
		Utah	EC-SM1-UT	650-203	EVDR-UT
		Medex Logical	EC-SM1-MX1		EVDR-ML
Infinity, Gamma X, XI Delta Delta XI Kappa, Kappa xlt Vista XI					Direct cables change above to EVDR2-XX (XX - transducer)

Drager IBP Adapter (Continued)


Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Drager Infinity Vista, Delta Kappa		Argon	EC-SM2-AG		EVDR1-1R
		B. Braun	EC-SM2-BB		EVDR1-BB
		BD	EC-SM2-BD	MS22148, 3375933	ERDR1-BD
		Edwards	EC-SM2-ED	MS22147	ERDR1-ED
		Medex Abbott	EC-SM2-MX	MS22535	ERDR1-DA
		Utah	EC-SM2-UT	MS22534	ERDR1-UT
		Medex Logical	EC-SM2-MX1		ERDR1-ML

EDAN (See 6-Pin IBP Adapter)


Marquette (GE) IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
GE Marquette Eagle, Dash 1000, Dash 2000, Dash 3000, Dash 3100, Dash 4000, Solar 7000, Solar 8000, Solar 9000, Tram Module 100-600, Cath Lab, MAC-Lab, Tram Scope 100-800 A or L Series, Tram Rac 3, Tram Rac 4A; Datex Ohmeda AM Module E-PRESTN		Argon	EC-MQ-AG		EVMA-AR
		B. Braun	EC-MQ-BB	22339301	EVMA-BB
		BD	EC-MQ-BD	700077-001, 684102	EVMA-BD
		Edwards	EC-MQ-ED	2021197-001	EVMA-ED
		Medex Abbott	EC-MQ-MX	700075-001	EVMA-MA
		Utah	EC-MQ-UT	700078-001, 650-236	EVMA-UT
		Medex Logical	EC-MQ-MX1		EVMA-ML

Hellige (GE) IBP Adapter


Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
SMK Series SMV Series SMU series		Argon	EC-HL-AG		EVHE-AR
		B. Braun	EC-HL-BB		EVHE-BB
		BD	EC-HL-BD		EVHE-BD
		Edwards	EC-HL-ED		EVHE-ED
		Medex Abbott	EC-HL-MX		EVHE-MA
		Utah	EC-HL-UT		EVHE-UT
		Medex Logical	EC-HL-MX1		EVHE-ML

Hwatime IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
H6, H8, H9 –see 6 pin HT6, HT8, HT9		Argon			EVHW-AR
		B. Braun			EVHW-BB
		BD			EVHW-BD
		Edwards			EVHW-ED
		Medex Abbott			EVHW-MA
		Utah			EVHW-UT
		Medex Logical			EVHW-ML

**Infinium, Invivo, Ivy,
Kontron, Medtronic,
MEK. Mennen**

Infinium IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Arrow KAAT I Plus and II Plus; DRE Waveline and Vida; Infinium Procyon and Vega; Kontron Minimon, Supermon, Fetalmon 2B Fetalogic 4000, K2000, IAB Balloon Pump, 7000 Series, Aero KAAT II Plus		Argon	EC-KTN-AG		EVIN-AR
		B. Braun	EC-KTN-BB		EVIN-BB
		BD	EC-KTN-BD		EVIN-BD
		Edwards	EC-KTN-ED		EVIN-ED
		Medex Abbott	EC-KTN-MX		EVIN-MA
		Utah	EC-KTN-UT	650-230	EVIN-UT
		Medex Logical	EC-KTN-MX1		EVIN-ML


Invivo (See 6-Pin IBP Adapter)

IVY Biomedical (See 6-Pin IBP Adapter)


Kontron (See Infinium)

Medtronic (See 6-Pin IBP Adapter)


MEK IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
		Argon	EC-MEK-AG		EVME-AR
		B. Braun	EC-MEK-BB		EVME-BB
		BD	EC-MEK-BD		EVME-BD
		Edwards	EC-MEK-ED		EVME-ED
		Medex Abbott	EC-MEK-MX		EVME-MA
		Utah	EC-MEK-UT		EVME-UT
		Medex Logical	EC-MEK-MX1		EVME-ML

Mennen IBP Adapter (Also see 6-Pin IBP Adapter)



Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Mennen Envoy and Vitalogik 5500		Argon	EC-MN-AG		EVMEN-AR
		B. Braun	EC-MN-BB		EVMED-BB
		BD	EC-MN-BD	684248	EVMEN-BD
		Edwards	EC-MN-ED		EVMEN-ED
		Medex Abbott	EC-MN-MX		EVMEN-MA
		Utah	EC-MN-UT		EVMEN-UT
		Medex Logical	EC-MN-MX1		EVMEN-ML

Mindray IBP Adapter (Also see 6-Pin IBP Adapter)


Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Beneview T5, T6, T8, DPM6, DPM7		Argon	EC-MR-AG		EVMI-AR
		B. Braun	EC-MR-BB		EVMI-BB
		BD	EC-MR-BD	001C-30-70757	EVMI-BD
		Edwards	EC-MR-ED	0010-21-12179	EVMI-ED
		Medex Abbott	EC-MR-MX		EVMI-MA
		Utah	EC-MR-UT	650-206	EVMI-UT
		Medex Logical	EC-MR-MX1		EVMI-ML

Mediana (See 6-Pin IBP Adapter)

Nihon Kohden IBP Adapter

Models	Adapter	Transducers	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Life Scope BSM: 1100 series, 1200 series, 1230A, 2100 series, 2300* series, 3100 series, 4100* series, 5100* series, 7100 series, 7200 series, 8300 series, 8500 series, 8800 series, 9500* series. Life Scope BSS 9800* series. Life Scope 8, 9, 11, 12, 14. Life Scope OEC 6000 and 7000 series. Life Scope OMP. *Requires adapter cable JP-910P from Nihon Kohden		Argon	EC-NK1-AG		EVNK-AR
		B. Braun	EC-NK1-BB		EVNK-BB
		BD	EC-NK1-BD	684090	EVNK-BD
		Edwards	EC-NK1-ED		EVNK-ED
		Medex Abbott	EC-NK1-MX		EVNK-MA
		Utah	EC-NK1-UT	650-225	EVNK-UT
		Medex Logical	EC-NK1-MX1		EVNK-ML
Life Scope TR, BSM-6000 series, BSM-6301K, BSM-6501K, BSM-6701K, BSM-9101K, Life Scope J		Argon	EC-NK2-AG		EVNK2-AR
		B. Braun	EC-NK2-BB		EVNK21-BB
		BD	EC-NK2-BD	JP-900P	EVNK2-BD
		Edwards	EC-NK2-ED	JP-902P	EVNK2-ED
		Medex Abbott	EC-NK2-MX		EVNK2-MA
		Utah	EC-NK2-UT		EVNK2-UT
		Medex Logical	EC-NK2-MX1		EVNK2-ML

Philips IBP Adapter


Models	Adapter	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
78341A, 78342A, 78352A, M1165A, CMS 24, V26, M2, M3, M4, M1026A, M1166A, M1175A, M1205A, M1350A, M1353A, 8040A, 3536A, all HP Monitors # XX, MXXXXA M1006A - BPM M1006B, Merlin 1300, Omnicare to M1176A		For Argon	EC-HP-AG		EVPH-AR
		For B. Braun	EC-HP-BB	M1634A	EVPH-BB
		For BD Transducers	EC-HP-BD	684081	EVPH-BD
		For Edwards	EC-HP-ED		EVPH-ED
		For Medex Abbott	EC-HP-MX		EVPH-MA
		For Utah Transducers	EC-HP-UT	650-206	EVPH-UT
		For Medex Logical	EC-HP-MX1		EVPH-ML

**Physio Control,
Siemens, Spacelabs,
Welch Allyn, Zoll,
Converter Cables**

Physio Control (See 6-Pin IBP Adapter)

Siemens (See Drager)

Spacelabs IBP Adapter (Also see 6-Pin IBP Adapter)





Models	Adapter	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
Ultraview		For Argon	EC-SL-AG		EVSP-AR
		For B. Braun	EC-SL-BB		EVSP-BB
		For BD Transducers	EC-SL-BD		EVSP-BD
		For Edwards	EC-SL-ED		EVSP-ED
		For Medex Abbott	EC-SL-MX		EVSP-MA
		For Utah Transducers	EC-SL-UT		EVSP-UT
		For Medex Logical	EC-SL-MX1		EVSP-ML

Welch Allyn (See 6-Pin IBP Adapter)




Zoll (See 6-Pin IBP Adapter)

CONVERTER CABLES

Drager IBP Converter



Specification	Converter	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
IBP converter cable used for converting a Spacelabs IBP connector to two Siemens 10-Pin adapter cables. 0.3m cable.		IBP Converter Cable	EC-S-SM1/2	5731281	EVSP-DRA
					
IBP converter cable used for converting a Spacelabs IBP connector to two Siemens 6-Pin adapter cables. 0.3m cable.		IBP Converter Cable	EC-S-SM2/2	5592147	EVSP-DR1A
					

Drager IBP Converter (Continued)



Specification	Converter	Options	Unimed/Evotec	ORIGINAL OEM #	EVOTEC OEM
IBP converter cable used for converting a Spacelabs IBP connector to two Siemens 10-Pin adapter cables. 0.3m cable.		IBP Converter Cable	E-SM2-SM1/2		EVSP-S11/2
IBP converter cable used for converting a Spacelabs IBP connector to a Siemens 6-Pin adapter cable. 0.3m cable.		IBP Converter Cable	E-SM1-SM2		EVSP-S12
IBP converter cable used for converting a Spacelabs IBP connector to a Siemens 10-Pin adapter cable. 0.3m cable.		IBP Converter Cable	E-SM2-SM1		EVSP-S11

TRANSDUCERS



Abbott Medex Disposable IBP Transducer

Photo	Connector	Options	Unimed/Evotec	EVOTEC OEM
		Box of 20 disposable IBP transducers	EP-MX	EVME-20



BD Disposable IBP Transducer

Photo	Connector	Options	Unimed/Evotec	EVOTEC OEM
		Box of 20 disposable IBP transducers	EP-BD	EVBD-20



Edwards TruWave Disposable IBP Transducer

Photo	Connector	Options	Unimed/Evotec	EVOTEC OEM
		Box of 20 disposable IBP transducers	EP-ED	EVED-20

B. Braun Disposable IBP Transducer

Photo	Connector	Options	Unimed/Evotec	EVOTEC OEM
		Box of 20 disposable IBP transducers	EP-BB	EVBB-20

Utah Disposable IBP Transducer

Photo	Connector	Options	Unimed/Evotec	EVOTEC OEM
		Box of 20 disposable IBP transducers	EP-UT-30	EVUT-20